

Original Research

The value of the SmartBx[™] system in improving the detection of prostate cancer in patients undergoing Transrectal Ultrasound-Guided Biopsy

Ismail Masarwa 1,† , Zaher Bahouth 1,†,* , Alexander Shefler 1 , Jacob Bejar 2 , Sagi Shprits 1 , Ofir Avitan 1 , Sarel Halachmi 1 and Ofer Nativ 1

DOI:10.31083/j.jmcm.2020.02.728

This is an open access article under the CC BY 4.0 license (https://creativecommons.org/licenses/by/4.0/).

Background: Prostate cancer is one of the most common cancers affecting men worldwide with an increasing incidence in most countries. The gold standard for diagnosis is Transrectal Ultrasonography-guided prostate biopsy. The quality of the biopsy is determined by parameters like the length of the biopsy core and the degree of fragmentation, affecting detection rate. SmartBx™ is a tool designed for the handling of prostate biopsy specimens and is intended for maximal tissue preservation. The aim of the current study was to evaluate the detection rate of prostate cancer using the new SmartBx™ tool. Patients and Methods: The study group had a biopsy using the SmartBx™ system and each biopsy core was applied into a cassette, which included two cores of biopsies which were then downloaded into a formalin-filled container. The control group included patients who underwent a Transrectal Ultrasonography-guided biopsy using the "standard" technique of biopsy collection. Data were collected prospectively in the study group, and retrospectively in the control group. Overall detection rate, and detection rate per core, were calculated and compared between the two groups. Results: The study included 600 patients, 300 in each group. There were no statistically significant differences between the groups regarding baseline characteristics. The average core numbers obtained per patient were 13.0 ± 2.7 and 14.1 ± 4.4 (P = 0.001) in the SmartBx[™] and control groups, respectively. The use of the SmartBx[™] system resulted in 21.2% better tissue preservation than in the control group. The overall detection rate as well as the detection rate per core were higher in the study group. More Gleason score 4 + 3 (Grade Group III) were detected in the study group (44.1% vs. 12.5%, P = 0.01). Conclusions: The use of the SmartBxTM system yielded more biopsy tissue for diagnosis. It increased the detection rate per core and the detection of Gleason score 4 + 3 (Grade Group III) prostate cancer.

Keywords

Prostate cancer; transrectal ultrasound; SmartBX

1. Introduction

Prostate cancer (PCa) is one of the most common cancers affecting men worldwide with an increasing incidence in most countries [1, 2]. Early diagnosis and treatment of PCa can reduce metastatic disease and mortality [3]. Major efforts have been invested to improve the detection rate of clinically significant PCa, and to eventually administer the treatment of choice [4, 5].

Transrectal ultrasound (TRUS)-guided prostate biopsy is the gold standard procedure for the pathological diagnosis of PCa [6, 7]. The quality of the biopsy is determined by parameters like the length of biopsy core and degree of fragmentation, affecting detection rate [8, 9, 10, 11]. It was also found that improved preembedding method by placing the biopsy core directly in cassettes may have some advantage if this is performed at the biopsy suite by the urologists [8]. This is mainly due to decreasing manipulation and damage to the biopsy specimen.

Poor biopsy quality may affect prognostic parameters like number of involved cores and extent of cancer in an individual core, which may affect treatment decision [12].

SmartBx $^{\text{TM}}$ (UC-Care Medical Systems, Yokneam, Israel) is a tool designed for the handling of prostate biopsy specimens and is intended for maximal tissue preservation as well as for accurate localization and orientation of the biopsy cores.

The main objective of the current study was to evaluate the impact of the SmartBx $^{\text{TM}}$ system on the detection rate of PCa.

2. Patients and methods

SmartBxTM is a semiautomatic device which enables downloading the biopsy core from the biopsy needle while maintaining its unity and orientation (Fig. 1). The SmartBxTM cassette comprises a Nylon-based membrane which allows non-covalent binding with the biopsy tissue via hydrophobic and charge-based interactions. This binding does not result in any modifications to the

¹Department of Urology, Bnai-Zion Medical Center, Haifa, Israel

²Department of Pathology, Bnai-Zion Medical center, Haifa, Israel

[†]These authors contributed equally.

^{*}Correspondence: Zaher.bahouth@b-zion.org.il (Zaher Bahouth)

tissue. The membrane carrying the biopsy cores is designed to go through the standard fixation and embedding processes.

The study was approved by the local ethical committee and all patients signed a written informed consent (Institutional Review Board no. BNZ-0053-11). We performed a prospective data collection for the study group which included patients with clinically suspected PCa that were referred to TRUS-guided biopsy in our hospital from 2011 to 2017. The results were compared to a control group who underwent TRUS-guided biopsy during the same years without the use of SmartBxTM, and their data were collected retrospectively.

In the study group, biopsies obtained by SmartBx TM were downloaded by the urologist into a formalin-filled container, with each cassette containing two cores. In the control group, biopsy cores were downloaded by shaking the needle directly into a formalin-filled container and two floating cores were preserved in each container.

All the biopsies were performed by a single urologist and sent to the local pathology lab and managed by the same standard tissue processing, and diagnosed by the same pathologists.

The detection rate was defined as the number of patients diagnosed with PCa divided by the total number of the patients in each group. Moreover, because of variation in the total number of cores, we calculated the detection rate per core by dividing the positive cores on the total number of the cores in each group.

The initial core length was documented on the needle notch using uEye camera (Imaging Development Systems, Obersulm, Germany), and measured again on the slide. From these measurements we calculated the percentage of biopsy tissue preservation. T-test was used to compare between continuous variables, and Chi-




Fig. 1. The SmartBxTM needle with the gun.

Square was used to compare between categorical variables. Statistically significant differences were defined as p value of < 0.05. All statistical analyses were performed using the SAS system.

3. Results

Six-hundred patients were included in the study, 300 in each group. Table 1 shows that there were no statistically significant differences in baseline characteristics between the two groups including age, PSA, digital rectal exam and prostate volume.

The average core number obtained per patient was 13.0 ± 2.7 and 14.1 ± 4.4 (P = 0.001) in the SmartBxTM and control groups, respectively. The use of SmartBxTM system resulted in 1.13 ± 4.07 mm longer tissue core (12.00 ± 2.80 mm vs. 10.87 ± 2.96 mm, P = 0.001) available for pathological evaluation compared with the standard system. This indicated a 21.2% better biopsy specimen preservation (84.5% vs. 63.3%).

The detection rate in the study group was higher than in the

Table 1. Baseline clinical characteristics of the entire cohort.

	$\mathrm{SmartBx^{TM}}$ $(n=300)$	${\rm Control}\ (n=300)$	P value
Mean age (years)	65.2 ± 8.5	65.8 ± 7.6	0.35
PSA (ng/mL)	8.8 ± 10.5	9.0 ± 9.5	0.8
PSA density (ng/mL ²)	0.17 ± 0.24	0.15 ± 0.17	0.23
Prostate volume (cc)	65.0 ± 36.0	68.8 ± 40.9	0.2
Abnormal DRE	45.3%	47.3%	0.6

PSA - Prostate-specific antigen. DRE -- Digital rectal examination.

Table 2. Pathological results of SmartBx $^{\text{TM}}$ group and control group.

	SmartB x^{TM} (n = 300)	Control ($n = 300$)	P value
Average core number per patient	13.0 ± 2.7	14.1 ± 4.4	0.001
Biopsy length on the slide (mm)	12 ± 2.8	10.87 ± 2.96	0.001
Percentage of biopsy preservation	84.5%	63.3%	0.01
Detection rate	37.7%	30.3%	0.06
Detection rate per core	(548/3892) 14.1%	(342/4239) 8.1%	0.01
Detection rate by GS (GG)			
6 (I)	33.0% (n = 37)	37.8 (n = 34)	0.48
7 = 3 + 4 (II)	21.4% (n = 24)	38.9% (n = 35)	0.01
7 = 4 + 3 (III)	17.0% (n = 19)	5.6% (n = 5)	0.01
8 (IV)	10.7% (n = 12)	19.0% (n = 9)	0.08
9-10 (V)	17.9% (n = 20)	7.8% (n = 7)	0.04

GS = Gleason Score. GG = ISUP Grade Group.

34 Masarwa et al.




Fig. 2. Two blopsy cores. The image on the right shows the needle crib. The left image shows two cores on a cassette.

control group (37.7% vs. 30.3%) although statistical significance was not attained (P = 0.06). The detection rate per core was substantially higher in the SmartBxTM group compared with the control group (14.1% vs. 8.1%, P = 0.001), as shown in Table 2. The latter Table also shows the differences in the detection rates stratified by Gleason score (GS) and Grade groups (GG). It should be noted that higher rates of GS 4 + 3 (GG III) were detected in the study group (44.1% vs. 12.5%, P = 0.01).

4. Discussion

Although being the gold-standard procedure for diagnosing PCa, TRUS-guided prostate biopsy has a detection rate of less than 50% [13, 14]. Every effort should be invested to maximize the detection rate, including proper processing of the specimen, which have been shown to improve its quality for proper pathological examination [15, 16].

Core length is proportionally related to better diagnosis and higher detection rate [9, 17]. Specimen observation is also related to improved pathological examination [8].

In our current study, we showed that using SmartBxTM system resulted in a 1.13 mm longer core in average, and 21.2% more biopsy tissue was preserved when compared to standard biopsy. This was eventually translated into higher detection rate per core. It should be mentioned that the overall detection rate, which was higher in the study group, did not reach statistically significance, though it was nearly attained (P = 0.06). This could be explained by the fact that the standard group had more cores than the study group.

One more interesting fact was the higher proportion of GS 4 + 3 (GG III) in the study group and the higher proportion of GS 3 + 4 (GG II) in the control group. We do not know the reason for this difference, however, it could be related to the fact that the longer and better-preserved tissue in the SmartBxTM group enabled the pathologists to detect a higher proportion of Gleason pattern 4, which differentiates between GS 3 + 4 (GG II) and 4 + 3 (GG III). This is important because of the different prognosis of these two groups, which differentiates patients between favorable and unfavorable intermediate risk PCa by NCCN tables [7].

Our study has some limitations, the main of which was the fact

that there was no randomization. Moreover, data on the control group was obtained retrospectively. The number of patients in each group is relatively small, and as such could miss other statistically significant differences between the groups. However, our study strength points were the facts that all the biopsies were performed by the same surgeon and all the specimens were examined by the same pathologist. Moreover, the data on the study group were collected prospectively.

In conclusion, this study showed that the use of SmartBxTM system preserved more tissue and core-length, improved detection rate per core, and increased detection rate of GS 4 + 3 (GG III).

Acknowledgments

We would like to express my gratitude to all those who helped me during the writing of this manuscript.

Conflict of interest

The authors declare no conflict of interest.

Submitted: May 01, 2020 Accepted: June 01, 2020 Published: June 20, 2020

References

- [1] Ferlay J, Soerjomataram I, Dikshit R, Eser S, Mathers C, Rebelo M, et al. Cancer incidence and mortality worldwide: sources, methods and major patterns in GLOBOCAN 2012. Int J Cancer. 2015; 136(5): E359–E386.
- [2] Center MM, Jemal A, Lortet-Tieulent J, Ward E, Ferlay J, Brawley O, et al. International variation in prostate cancer incidence and mortality rates. Eur Urol. 2012; 61(6): 1079–1092.
- [3] Schröder FH, Hugosson J, Carlsson S, Tammela T, Määttänen L, Auvinen A, et al. Screening for prostate cancer decreases the risk of developing metastatic disease: findings from the European Randomized Study of Screening for Prostate Cancer (ERSPC). Eur Urol. 2012; 62(5): 745–752.
- [4] Zhou Steve R, Edward C, Aalhad P, Jiaoti H, Marks Leonard S, Shyam N. Prostate cancer detection rate of freehand versus 3-Dimensional template mapping biopsy using a magnetic resonance imaging-ultrasound fusion device in biopsy Naïve men. *J Urol.* 2020; 203(4): 699-705.
- [5] Campa R, Del Monte M, Barchetti G, Pecoraro M, Salvo V, Ceravolo I, et al. Improvement of prostate cancer detection combining a computer-aided diagnostic system with TRUS-MRI targeted biopsy. *Abdom Radiol.* 2019; 44(1): 264–271.
- [6] Mottet N, Bellmunt J, Bolla M, Briers E, Cumberbatch MG, De Santis M, et al. EAU-ESTRO-SIOG guidelines on prostate cancer. Part 1: screening, diagnosis, and local treatment with curative intent. Eur Urol. 2017; 71(4): 618–629.
- [7] National Comprehensive Cancer Network. Prostate cancer (Version 4.2019) [Internet]. [cited 2020 Feb 1]. Available from: https://www.nccn.org/professionals/physician gls/pdf/prostate.pdf.
- [8] Van der Kwast T, Bubendorf L, Mazerolles C, Raspollini MR, Van Leenders GJ, Pihl CG, et al. Guidelines on processing and reporting of prostate biopsies: the 2013 update of the pathology committee of the European Randomized Study of Screening for Prostate Cancer (ERSPC). Virchows Arch. 2013; 463(3): 367–377.
- [9] Iczkowski KA, Casella G, Seppala RJ, Jones GL, Mishler BA, Qian J, et al. Needle core length in sextant biopsy influences prostate cancer detection rate. *Urology*. 2002; 59(5): 698–703.
- [10] Doluoglu OG, Yuceturk CN, Eroglu M, Ozgur BC, Demirbas A, Karakan T, et al. Core length: an alternative method for increasing cancer detection rate in patients with prostate cancer. Urology journal. 2015; 12(5): 2324–2328.
- [11] Ergun M, Islamoglu E, Yalcinkaya S, Tokgoz H, Savas M. Does

Volume 3, Number 2, 2020 35

- length of prostate biopsy cores have an impact on diagnosis of prostate cancer? *Türk Ürol DergisiTurkish J Urol*. 2016; 42(3): 130–133
- [12] Steyerberg EW, Roobol MJ, Kattan MW, van der Kwast TH, de Koning HJ, Schröder FH. Prediction of Indolent Prostate Cancer: Validation and Updating of a Prognostic Nomogram. J Urol. 2007; 177(1): 107–112
- [13] Bhindi B, Mamdani M, Kulkarni GS, Finelli A, Hamilton RJ, Trachtenberg J, et al. Impact of the U.S. Preventive Services Task Force recommendations against prostate specific antigen screening on prostate biopsy and cancer detection rates. J Urol. 2015; 193(5): 1519–1524.
- [14] Ravery V, Billebaud T, Toublanc M, Boccon-Gibod L, Hermieu JF, Moulinier F, et al. Diagnostic value of ten systematic TRUS-Guided

- prostate biopsies. Eur Urol. 1999; 35(4): 298-303.
- Boccon-Gibod L, van der Kwast ThH, Montironi R, Boccon-Gibod L, Bono A. Handling and pathology reporting of prostate biopsies. *Eur Urol.* 2004; 46(2): 177–181.
- [16] Tombal B, Tajeddine N, Cosyns J-P, Feyaerts A, Opsomer R, Wese FX, et al. Does site-specific labelling and individual processing of sextant biopsies improve the accuracy of prostate biopsy in predicting pathological stage in patients with T1c prostate cancer?: SITE-SPECIFIC LABELLING OF BIOPSY IN PROSTATE CANCER. BJU Int. 2002; 89(6): 543–548.
- Obek C, Doğanca T, Erdal S, Erdoğan S, Durak H. Core length in prostate biopsy: size matters. J Urol. 2012; 187(6): 2051–2055.

36 Masarwa et al.